



704.CELLULAR IMMUNOTHERAPIES: EARLY PHASE AND INVESTIGATIONAL THERAPIES |
NOVEMBER 28, 2023

Varnimcabtagene Autoleucel (IMN-003A): Differences in T Cell Subset Phenotype in B-ALL and B-NHL Cohorts Did Not Influence Efficacy Outcomes in the Phase-2 Study (IMAGINE)

Sunil Bhat, Sharat Damodar, Raja Thirumalairaj, Pankaj Malhotra, Alex Jose Pulikkottil, Jeetendra Kumar, Melina Soares, Sudeshna Dhar, Sandeepa Shetty, Mani Bharathi T.I., Gopinadh Jakka, Anne Roshan Joseph, Arun Kumar MG, Pallavi Arasu, Sri Ramulu Elluru, Mohammed Manzoor Akheel, Giridhar K Chenji, Lakshmikanth Gandikota, Rahul Nahar, Anil Kamat



Check for updates

Blood (2023) 142 (Supplement 1): 6881.

<https://doi.org/10.1182/blood-2023-184521>



Split-Screen



Share 



Tools 

Background: Varnimcabtagene autoleucel (IMN-003A) is an autologous CD19 directed CAR-T cell product with a 4-1BB co-stimulatory domain and a non-FMC63 murine single chain variable fragment (A3B1 binder), manufactured in India, tested in the IMAGINE study (CTRI/2022/03/041162), a phase-2 clinical trial for patients (pts) with relapsed and/or refractory B cell malignancies. A fractionated infusion of 1×10^6 CAR+ cells (IMN-003A)/kg for B-ALL cohort and 5×10^6 CAR+ cells (IMN-003A)/kg for B-NHL cohort was administered over 3 days (10%, 30%, 60%). Phenotypic data comparing the apheresis (AP) and final product (FP) from B-ALL and B-NHL cohorts, related manufacturing and clinical response observations are discussed in this abstract.

Methods: IMN-003A was manufactured using the Miltenyi CliniMACS Prodigy[®], a cGMP compliant closed system. The apheresis and final infusion product (FP) were analysed for percent CAR transduction and T cell subsets by flow cytometry. CAR copies in the final product were determined by ddPCR. Cell counts measured using automated Vi-cell counter were used to calculate fold expansion.

Results: Apheresis PBMC from 25 patients in the IMAGINE study were selected for T-cells, transduced & expanded to achieve the FP, which was infused in 24 patients (1 patient withdrawn).

[Skip to Main Content](#)

Three patients needed a second apheresis and re-manufacturing to achieve the target dose. The FP comprised primarily of CD3+ T cells (median 99.38% range 96.3 - 99.92%) with mean of



CAR+ T-cells for B-NHL cohort (range: 8.5 - 51.07%, $p=0.02$). The vector copy number (VCN) however did not differ between the two cohorts (mean of 2.66 and 2.37 copies for B-ALL and B-NHL cohorts respectively; $p=0.4$). The fold expansion at Day 6 of manufacturing was also significantly lower for B-NHL cohort (6.93) compared to the B-ALL cohort (10.33, $p=0.02$). Consistent with the relatively more challenging manufacturing for the B-NHL cohort with 5 times higher target dose, it was observed that the mean proportion of naïve T-cells ($CCR7^+ CD45RA^+$) was lower in both AP (B-NHL 9.71%, B-ALL 27.24%, $p<0.01$) and FP (B-NHL 16.68%, B-ALL 48.94%, $p<0.0001$), while the terminally differentiated T_{EMRA} and effector memory T_{EM} subsets were higher in the AP (B-NHL 53.8%, B-ALL 45.08%, $p=0.175$) & FP (B-NHL 48.77%, B-ALL 29.54%, $p=0.021$), highlighting the T-cell phenotypic differences in the two cohorts.

The mean CD4/CD8 ratio in AP & FP was 0.56 and 1.14 respectively. There was reversal of CD4/CD8 ratio from apheresis to final product (CAR+ T cells, ratio=3.37) which was statistically significant ($p<0.01$). The CD4:CD8 ratio did not differ between B-ALL and B-NHL cohorts. This could potentially be explained by the higher proportion of naïve T-cells observed within CD4 subset compared to CD8 subset as against the proportion of T_{EMRA} cells. Despite the higher fraction of CD4 CAR-T cells within the final product infused, post-infusion, there was a statistically significant trend ($p<0.01$) towards reversal of CD4/CD8 ratio in the peripheral blood from D0 to D+90. Ongoing flow-cytometry based analysis shall help further characterise the CAR+ T-cells post-infusion.

Interestingly, despite phenotypic differences in T-cell subsets between B-ALL and B-NHL, overall response rate (ORR) at D+90 was very similar for both cohorts [B-ALL 80% ($n=8/10$); B-NHL 81.8% (9/11)]. Updated results shall be presented at the meeting.

Conclusions: Differences in the manufacturing outcomes (CAR % transduction, fold expansion) and T cell phenotype (naïve and effector T_{EMRA} cell subsets) were observed in Apheresis and FP, between the B-ALL and B-NHL cohorts. A higher fraction of relatively more differentiated T cell phenotypes of effector (T_{EMRA}) and effector memory (T_{EM}) and a lower fraction of naïve T cells in the B-NHL cohort and within CD8 T cells, reflects a manufacturing conundrum for B-NHL patients with a higher target dose and the CD4:CD8 ratio in the final product. However, these differences in the final product did not adversely impact the clinical responses at D+90 with equivalent efficacy in both cohorts. The ideal T cell phenotype or their proportions within the final product for maximum efficacy thus remains uncertain. A deeper characterisation including exhaustion status of the apheresis, final product and post-infusion samples for the T cell subsets can provide further insights.

Disclosures

Pulikkottil: Immuneel Therapeutics Private Limited: Current Employment. **Kumar:** Immuneel Therapeutics Private Limited: Current Employment. **Soares:** Immuneel Therapeutics Private Limited: Current Employment. **Dhar:** Immuneel Therapeutics Private Limited: Current Employment.

Employment. **Joseph:**Immuneel Therapeutics Private Limited: Current Employment. **Kumar**  **MG:**Immuneel Therapeutics Private Limited: Current Employment. **Arasu:**Immuneel Therapeutics Private Limited: Current Employment. **Elluru:**Immuneel Therapeutics Private Limited: Current Employment. **Akheel:**Immuneel Therapeutics Private Limited: Current Employment. **Chenji:**Immuneel Therapeutics Private Limited: Current Employment. **Gandikota:**Immuneel Therapeutics Private Limited: Current Employment. **Nahar:**Immuneel Therapeutics Private Limited: Current Employment. **Kamat:**Immuneel Therapeutics Private Limited: Current Employment.

Figure 1




Figure 1. T cell phenotypic differences in apheresis and FP between B-ALL and B-NHL patients: (A) Representative flow cytometry plots from two patients each from B-ALL & B-NHL cohorts, showing the different populations (Naïve T cells, T_{CM}, T_{EM} and T_{EMRA}) for Apheresis and the Final Product. Parent gate is live CD45+CD3+ cells. (B) Dot plots showing the differences in Naïve T cells and T_{EMRA} phenotype between B-ALL and B-NHL cohorts for both Apheresis and the Final Product.

[VIEW LARGE](#) [DOWNLOAD SLIDE](#)

[Sign in via your Institution](#)

[Add comment](#)

[Skip to Main Content](#)

